Quantibody[®] Human Ig Isotype Array 1

--Quantitative measurement of 8 human Immnunoglobulins

Patent Pending Technology

User Manual (Version July 2010)

Cat # QAH-ISO-1



We Provide You With Excellent Protein Array Systems and Service

Tel:(Toll Free) 1-888-494-8555 or 770-729-2992; Fax: 1-888-547-0580; Website:<u>www.raybiotech.com</u> Email: <u>info@raybiotech.com</u>

Immunoglobulin Detected (8)	IgA, IgD, IgE, IgM, IgG1, IgG2, IgG3, and IgG4
Format	One standard glass slide is spotted with 16 wells of identical Immunoglobulin antibody arrays. Each antibody is arrayed in quadruplicate.
Detection Method	Fluorescence with laser scanner: Cy3 equivalent dye
Sample Volume	$50 - 100 \ \mu l \ per \ array$
Serum Dilution	40,000x - 200,000x
Reproducibility	CV <20%
Assay duration	6 hrs



Quantibody® Human Immunoglobulin Array 1 1

TABLE OF CONTENTS

I.	Overview	1
	Introduction	3
	How It Works	6
II.	Materials Provided	7
	Additional Materials Required	7
III.	General Considerations	8
	A. Preparation of Samples	8
	B. Handling Glass Chips	8
	C. Incubation	8
IV.	Protocol	9
	A. Complete Air Dry the Glass Chip	9
	B. Prepare Immunoglobulin Standard Dilutions	.9
	C. Blocking and Incubation	10
	D. Incubation with Detection Antibody Cocktail	11
	E. Incubation with Cy3 Equivalent Dye-Streptavidin	11
	F. Fluorescence Detection	11
	G. Data Analysis	12
V.	Immunoglobulin Array Map & Standard Curves	13
VI.	8-Point Standards	14
VII.	System Flow	15
VIII.	Quantibody [®] Q-Analyzer	16
IX.	Troubleshooting Guide	17
Х.	Select Quantibody [®] Publications	18
XI.	Experimental Record Form	19
XII.	How to Choose Quantibody [®] Products	20

I. Introduction

The human immune system consists of two functional components classified as the innate system (the physical, biochemical and cellular barriers), and the adaptive immune system (including lymphocytes and immunoglobulins). Immunoglobulins are the key elements of the humoral immune response in against polypeptide vertebrate parasitic invasion. The chains of immunoglobulins composed of two identical heavy (H) chains and two identical light (L) chains linked together by inter-chain disulfide bonds. While the amino-terminal portions that exhibits highly variable amino-acid composition are involved in antigen binding, the C terminal constant parts are involved in complement binding, placental passage and binding to cell membranes. Based upon the variation of the constant region of the heavy chain, nine immunoglobulin heavy chain isotypes are found in humans: IgA (with subclasses IgA1 and IgA2), IgD, IgE, IgM, and IgG (with subclasses IgG1, IgG2, IgG3, and IgG4).

IgG is the predominant immunoglobulin in the serum (about 12 mg/ml), which accounts for 75% of the total serum antibody of healthy individuals. IgG has a molecular weight of about 150 kDa. Four distinct subgroups of human IgG (IgG1, IgG2, IgG3, and IgG4) were first demonstrated in the 1960's by using polyclonal antisera prepared in animals immunized with human myeloma proteins. Quantitatively, the relative abundance of the four subclasses in adult human serum follows IgG1 > IgG2 > IgG3 = IgG4, which accounts for 6.98, 3.8, 0.56, and 0.56 mg/ml respectively.

IgA exists as a 160kd monomer in serum and as a 400kd dimer in secretions. Quantitatively, IgA is synthesized in amounts greater than IgG. However, due to its short half life in serum (6 days of IgA vs 21 days of IgG) and the lost of secretory form, the normal IgA serum level (2-3 mg/ml) ranked second after IgG, which accounts for 15% of the total antibody. There are two subclasses based on different heavy chains, IgA1 and IgA2. IgA1 is produced in bone marrow and makes up over 90% of the serum IgA. Secretory IgA is the predominant immunoglobulin present in gastrointestinal fluids, nasal secretions, saliva, tears and other mucous secretions of the body.

IgM is the third most common serum immunoglobulin (about 1.5 mg/ml) which makes up about 10% serum antibody. IgM normally exists as a pentamer (about 900 kDa) and has a theoretical valence of 10. As a consequence of its pentameric structure, IgM is a good antigen agglutinating and complement fixing immunoglobulin.

IgE exists as a 190 kDa monomer and is the least common serum immunoglobulins which accounts for 0.002% of the total serum antibodies. IgE is involved in allergic reactions. If an infectious agent succeeds in penetrating the IgA barrier, it comes up against the next line of defense, the IgE manned MALT (mucosa-associated lymphoid tissues) system. Contact with the allergen leads to the release of various pharmacological mediators that result in allergic symptoms.

IgD is a 175kd molecule that resembles IgG in its monomeric form. IgD is found in low level in serum (0.03 mg/ml) with uncertain serological functions. IgD antibodies are found for the most part on the surfaces of B lymphocytes. It is thought that IgD and IgM function as mutually-interacting antigen receptors for control of B-cell activation and suppression. Hence, IgD may have an immunoregulatory function.

The levels of different immunoglobulin subclasses follow a typical pattern in a healthy ethnic adult and are normally within a certain percentile ranges. Upon different antigenic stimulation, an antibody response will behave differently in the distribution of the different subclasses in plasma, such as increase, diminish or even the deficiency of producing one of the specific immunoglobulin subclass. Over the last decades numerous reports have appeared on the distribution of the immunoglobulin subclasses produced during immune responses to bacterial, viral, and parasitic antigens; autoantigens; tumor antigens, and many parenterally administered substances such as hormones, drugs, and allergens. As a result, quantification of the immunoglobulin isotype level in a given serum sample will provide the useful information about the myeloma states as well as in monitoring immunoglobulin replacement, plasmaphoresis, intravenous and immunosuppression therapy.

Quantitative measurement of the immunoglobulin subclasses can be done with Radial Immunodiffusion assay (RID), Nephelometry and turbidimetry assay, Radio Immuno Assay (RIA), Immuno-affnity chromatography, Direct Antiglobulin Test (DAT), or Enzyme-linked Immunosorbent Assay (ELISA). While most assays can detect only one subclass of the immunoglobulin a time, taking advantage of the array technology and the availability of the isotype specific monoclonal antibodies, Raybiotech Inc is proudly offering the research community with the Quantibody® Human Ig Isotype kit which can simultaneously and quantitatively detect multiple immunoglobulin subclasses in one experiment.

Quantibody[®] Human Ig Isotype Array uses sandwich-ELISA based technology for quantitative measurement of the eight human isotype immunoglobulins (IgG1, IgG2, IgG3, IgG4, IgA, IgD, IgE, and IgM) in human serum/plasma. Similar technology has been successfully used in our other quantibody[®] products for quantitative measurement of up to 40 cytokines in human, mouse, rat, and porcine samples. (See Section XI). Briefly, the 8 human immunoglobulin subclass-specific antibodies are arrayed in quadruplicate (together with two positive controls) with 16 identical sub-arrays in one standard glass slide. The kit also provides a myeloma-derived standard mixture of these 8 immunoglobulins, whose concentration has been predetermined.

In the experiment, standard immunoglobulins and samples are assayed in each well simultaneously through a sandwich like ELISA procedure. The signals will be detected using fluorescence-based detection method for consistency and reliability. By comparing signals from unknown samples to the standard curve generated for each of the 8 immunoglobulins, the unknown immunoglobulin concentration in the samples will be determined.

The kit provides a highly sensitive approach (within nano gram range) to simultaneously detect 8 immunoglobulin subclasses expression levels. The experimental procedure is simple and can be performed in any laboratory.

How It Works



I. Materials Provided

Upon receipt, all components of the Quantibody[®] Array kit should be stored at -20° C. At -20° C the kit will retain complete activity for up to 6 months. Once thawed, the glass chip, Immunoglobulin standard mix, detection antibody cocktail and Cy3 equivalent dye-conjugated Streptavidin should be kept at -20° C and all other components may be stored at 4° C. The entire kit should be used within 6 months of purchase.

Components

Item	Description	1-Slide kit	2-Slide kit
1	Quantibody [®] Array Glass Chip	1	2
2	Sample Diluent	1	1
3	20X Wash Buffer I	2	3
4	20X Wash Buffer II	1	1
5	Lyophilized Immunoglobulin standard mix *	1	1
6	Detection antibody cocktail	1	2
7	Cy3 equivalent dye-conjugated Streptavidin	1	2
8	Slide Washer/Dryer	1	1
9	Adhesive device sealer	5	10
10	Manual	1	1

* See Section VI for detailed Immunoglobulin concentrations after reconstitution.

Additional Materials Required

- Orbital shaker
- Laser scanner for fluorescence detection
- Aluminum foil
- Distilled water
- 1.5ml Polypropylene microcentrifuge tubes

II. General Considerations

A. Preparation of Samples

- Blood samples should be collected by venipuncture. Allow to clot naturally. Undiluted samples may be stored at 2-8 ^oC for up to 72 hours or in -20^oC for longer periods. Avoid repeated freezing and thawing.
- **Sample dilution**: The suggested dilution for the patient sample is 1:40,000. However, user may decide to use the optimum range for his own sample. Dilute 1ul serum sample in 199ul sample diluent. Gently mix well, and then proceed with another 1:200 dilution by adding 1ul of the diluted sample to 199ul sample diluent. The net dilution is 1:40,000.

B. Handling glass chips

- Do not touch the surface of the slides, as the microarray slides are very sensitive. Hold the slides by the edges only.
- Handle all buffers and slides with latex free gloves.
- Handle glass chip in clean environment.
- Because there is no barcode on the slide, transcribe the slide serial number from the slide bag to the back of the slide with a permanent marker before discarding the slide bag. Once the slide is disassembled, you might not have enough info to distinguish one slide from the other.

C. Incubation

- Completely cover array area with sample or buffer during incubation.
- Avoid foaming during incubation steps.
- Perform all incubation and wash steps under gentle rotation.
- Cover the incubation chamber with adhesive film during incubation, particularly when incubation is more than 2 hours or $<70 \ \mu$ l of sample or reagent is used.
- Several incubation steps such as step 6 (blocking), step 7 (sample incubation), step 10 (detection antibody incubation), or step 13 (Cy3 equivalent dye-streptavidin incubation) may be done overnight at 4^oC. Please make sure to cover the incubation chamber tightly to prevent evaporation.

III. Protocol

A. Completely air dry the glass chip

1. Take out the glass chip from the box, and let it equilibrate to room temperature inside the sealed plastic bag for 20-30 minutes. Remove slide from the plastic bag; peel off the cover film, and let it air dry at room temperature for another 1-2 hours.

Note: Incomplete drying of slides before use may cause the formation of "comet tails".

B. Prepare Immunoglobulin Standard Dilutions





- 2. Reconstitute the Immunoglobulin Standard Mix (lyophilized) by adding 500µl Sample Diluent to the tube. For best recovery, always quick-spin vial prior to opening. Dissolve the powder thoroughly by a gentle mix. Labeled the tube as Std1.
- 3. Label 6 clean microcentrifuge tubes as Std2 to Std7. Add 200µl Sample Diluent to each of the tubes.
- 4. Pipette 100µl Std1 into tube Std2 and mix gently. Perform 5 more serial dilutions by adding 100ul Std2 to tube Std3 and so on.

5. Add 100µl Sample Diluent to another tube labeled as CNTRL. Do not add standard Immunoglobulin or samples to the CNTRL tube, which will be used as negative control. For best results, include a set of standards in each slide.

Note: Since the starting concentration of each Immunoglobulin is different, their serial concentrations from Std1 to Std7 are varied which can be found in section VI.

C. Blocking and Incubation

- 6. Add 100µl Sample Diluent into each well and incubate at room temperature for 30 min to block slides.
- Decant buffer from each well. Add 100µl standard Immunoglobulins or samples to each well. Incubate arrays at room temperature for 1-2 hour. (*Longer incubation time is preferable for higher signals*)

Note: The sample volume can be 50-100 μ l. If sample volume is less than 70 μ l, cover the gasket with adhesive sealer to prevent evaporation during incubation. This step may be done overnight at 4^oC for best results.

Note: We recommend using a 1:40,000 dilution of the patients' serum samples. (See Section II for sample preparation.)

- 8. Wash:
 - Decant the samples from each well, and wash 5 times (5 min each) with 150 μ l of 1x Wash Buffer I at room temperature with gentle shaking. Completely remove wash buffer in each wash step. Dilute 20x Wash Buffer I with H₂O.
 - (*Optional*) Put the glass chip with frame into a box with 1x Wash Buffer I (cover the whole glass slide and frame with Wash Buffer I), and wash at room temperature with gentle shaking for 20 min.
 - Decant the 1x Wash Buffer I from each well, wash 2 times (5 min each) with 150 µl of 1x Wash Buffer II at room temperature with

Note: Incomplete removal of the wash buffer in each wash step may cause "dark spots". (Background signal is higher than that of the spot.)

D. Incubation with detection antibody cocktail and wash.

- 9. Reconstitute the detection antibody by adding 1.4 ml of Sample Diluent to the tube. Spin briefly.
- 10. Add 80 μl of the detection antibody cocktail to each well. Incubate at room temperature for 1-2 hour. (*Longer incubation time is preferable for higher signals and backgrounds*)
- 11. Decant the samples from each well, and wash 5 times with 150 μl of 1x Wash Buffer I and then 2 times with 150 μl of 1x Wash Buffer II at room temperature with gentle shaking. Completely remove wash buffer in each wash step.

E. Incubation with Cy3 equivalent dye -Streptavidin and wash

- 12. After briefly spinning down, add 1.4 ml of Sample Diluent to Cy3 equivalent dye-conjugated streptavidin tube. Mix gently.
- 13. Add 80 μl of Cy3 equivalent dye-conjugated streptavidin to each well. Cover the device with aluminum foil to avoid exposure to light or incubate in dark room. Incubate at room temperature for 1 hour.
- 14. Decant the samples from each well, and wash 5 times with 150 μl of 1x Wash Buffer I at room temperature with gentle shaking. Completely remove wash buffer in each wash step.

F. <u>Fluorescence Detection</u>

15. Disassemble the device by pushing clips outward from the slide side. Carefully remove the slide from the gasket. (Be careful not to touch the surface of the array side)



- 16. Place the slide in the slide Washer/Dryer (a 4-slide holder/centrifuge tube), add enough 1x Wash Buffer I (about 30 ml) to cover the whole slide, and then gently shake at room temperature for 15 minutes. Decant Wash Buffer I. Wash with 1x Wash Buffer II (about 30 ml) with gentle, and gently shake at room temperature for 5 minutes.
- 17. Remove water droplets completely by one of the following ways:
 - Put the glass chip into the Slide Washer/Dryer, and dry the glass chip by centrifuge at 1,000 rpm for 3 minutes without cap.
 - Or, dry the glass chip by a compressed N₂ stream.
 - Or gently apply suction with a pipette to remove water droplets. Do not touch the array, only the sides.
- 18. Imaging: The signals can be visualized through use of a laser scanner equipped with a Cy3 wavelength such as Axon GenePix. Make sure that the signal from the well containing the highest standard concentration (Std1) receives the highest possible reading, yet remains unsaturated.

Note: In case the signal intensity for different Immunoglobulin varies greatly in the same array, we recommend using multiple scans, with a higher PMT for low signal Immunoglobulins, and a low PMT for high signal Immunoglobulins.

G. Data Analysis

19. Data extraction can be done with most of the microarray analysis software (GenePix, ScanArray Express, ArrayVision, or MicroVigene). For quantitative data analysis, our Quantibody[®] Q-Analyzer software is available. It gives visual output as well as digital values. More information can be found in section VIII.

IV. Immunoglobulin Array Map & Standard Curves

POS1	POS2
IgA	IgD
IgE	IgM
IgG1	IgG2
IgG3	IgG4



V. 8-Point Standards

After reconstitution of the lyophilized Immunoglobulin standard mix, the 8point Immunoglobulin concentration used for generating the standard curve of a given antigen is listed below. The detection sensitivity of each protein in one experiment is user dependent. Try our array specific Quantibody Q-Analyzer to see your Limit of Detection (LOD). (Section VIII).

(pg/ml)	Cntrl	Std7	Std6	Std5	Std4	Std3	Std2	Std1
40 /								
IgA	0	274	823	2,469	7,407	22,222	66,667	200,000
	0	(0)	200	(17	1.052		1(((7	50.000
IgD	0	69	206	61/	1,852	5,556	16,667	50,000
IgE	0	69	206	617	1,852	5,556	16,667	50,000
IgM	0	137	412	1,235	3,704	11,111	33,333	100,000
IgG1	0	549	1,646	4,938	14,815	44,444	133,333	400,000
IgG2	0	549	1,646	4,938	14,815	44,444	133,333	400,000
IgG3	0	69	206	617	1,852	5,556	16,667	50,000
IgG4	0	69	206	617	1,852	5,556	16,667	50,000

Serial standard concentration (pg/ml)



VII. Quantibody® Q-Analyzer

Quantibody Q-Analyzer is an array specific, Excel-based program. However, it is not a simple calculation macro as it contains sophisticated data analysis.

Key features:

- <u>Simplicity:</u> Easy to operate and requires no professional training. With a simple copy and paste process, the cytokine concentration is determined.
- <u>Outlier Marking & Removing</u>: The software can automatically mark and remove the outlier spots for more accurate data analysis
- *<u>Normalization</u>*: The program allows for intra- and inter-slide normalization for large number of samples.
- <u>*Two Positive Controls*</u>: The program takes the two positive controls in each array for normalization.
- <u>*Two Analytical Algorithms*</u>: Users can choose either linear regression or log-log algorithms to meet their analytical needs.
- <u>*Two Data Outputs*</u>: standard curves and digital concentration.
- <u>User Intervention</u>: The program allows for user manual handling of those outliers and other analytical data.
- <u>Lower and Upper Limits Determination</u>: The program automatically marks out the values below or above the detection range.
- <u>Standard Deviation</u>: The program outputs the standard deviations of the quadruplicate spots for data accuracy.
- <u>Analytical Tips:</u> Q-Analyzer analysis tips are included in the program.

VIII. Troubleshooting guide

Problem	Cause	Recommendation
	Inadequate detection	Increase laser power and PMT parameters
	Inadequate reagent volumes or	Check pipettes and ensure correct
	improper dilution	preparation
	Short incubation time	Ensure sufficient incubation time and
Weak Signal		change sample incubation step to overnight
	Too low protein concentration in	Don't make too low dilution or concentrate
	sample	sample
	Improper storage of kit	Store kit as suggested temperature. Don't freeze/thaw the slide.
	Bubble formed during incubation	Avoid bubble formation during incubation
Uneven signal	Arrays are not completed covered by reagent	Completely cover arrays with solution
	Reagent evaporation	Cover the incubation chamber with adhesive
		film during incubation
	Cross-contamination from	Avoid overflowing wash buffer
	Comet tail formation	Air dry the slide for at least 1 hour before
		usage
	Inadequate standard reconstitution or	Reconstitute the lyophilized standard well at
	Improper dilution	the room temperature before making serial
Poor standard curve		dilutions. Check pipettes and ensure proper serial dilutions.
	Inadequate detection	Increase laser power that the highest
		standard concentration for each cytokine
		receives the highest possible reading yet
	Use freeze thewad exterine standards	remains unsaturated.
	Use neeze-mawed cytokine standards	new set of experiment Discard any leftover
	Overexposure	Lower the laser power
	Dark spots	Completely remove wash buffer in each
Iliah	-	wash step.
High	Insufficient wash	Increase wash time and use more wash
vackgi vullu		buffer
	Dust	Work in clean environment
	Slide is allowed to dry out	Don't dry out slides during experiment.

IX. Select Quantibody Publications

- Stechova, et al. Influence of Maternal Hyperglycaemia on Cord Blood Mononuclear Cells in Response to Diabetes-associated Autoantigens. *Scandinavian Journal of Immunology*. 2009. 70(2):149-158
- 2. Willingham, SB et al. NLRP3 (NALP3, Cryopyrin) facilitates in vivo caspase-1 activation, necrosis, and HMGB1 release via inflammasome-dependent and independent pathways. *J Immunol.* 2009; 183(3):2008-15
- 3. El Karim et al. Neuropeptides Regulate Expression of Angiogenic Growth Factors in Human Dental Pulp Fibroblasts. *Journal of Endodontics*, 2009; 35(6): 829-833
- Souquière S. et al. T-Cell tropism of simian T-cell leukaemia virus type 1 and cytokine profiles in relation to proviral load and immunological changes during chronic infection of naturally infected mandrills (*Mandrillus sphinx*). J Med Primatol. 2009; 38(4):279-89
- 5. Sharma, et al. Induction of multiple pro-inflammatory cytokines by respiratory viruses and reversal by standardized *Echinacea*, a potent antiviral herbal extract. *Antiviral Research*. 2009; 83(2)165-170.
- 6. Altamirano-Dimas, et al. *Echinacea* and anti-inflammatory cytokine responses: Results of a gene and protein array analysis. *Pharmacuetical Biology*. 2009; 47(6): 500-508.
- 7. Cheung, et al. Cordysinocan, a polysaccharide isolated from cultured *Cordyceps*, activates immune responses in cultured T-lymphocytes and macrophages: Signaling cascade and induction of cytokines. *Journal of Ethonopharmacology*. 2009; 124(1): 61-68.
- 8. Du, et al. P2-380: Identification and characterization of human autoantibodies that may be used for the treatment of prion diseases. *Alzheimer's and Dementia*. 2009; 4(4): T484-T484.
- 9. Van Rossum et al. Granulocytosis and thrombocytosis in renal cell carcinoma: a proinflammatory cytokine response originating in the tumour. *Neth J Med.* 2009; 67(5):191-4.
- Zhai, et al. Coordinated Changes in mRNA Turnover, Translation, and RNA Processing Bodies in Bronchial Epithelial Cells following Inflammatory Stimulation. *Molecular and Cellular Biology*. 2008; 28(24): 7414-7426.
- 11. Gao, et al. A Chinese herbal decoction, Danggui Buxue Tang, activates extracellular signal-regulated kinase in cultured T-lymphocytes. *FEBS Letters*, 2007; 581(26): 5087-5093. (This reference validates mulitplex ELISA results for several analytes with standard ELISA test results).
- 12. Piganelli, et al: Autoreactive T-cell responses: new technology in pursuit of an old nemesis. (Editorial Review) *Pediatric Diabetes* 2007: 8: 249–251

X. Experiment Record Form

Date: _____

File Name:

Laser Power: _____

PMT: _____

Well No.	Sample Name	Dilution factor
1	CNTRL	
2	Std7	
3	Std6	
4	Std5	
5	Std4	
6	Std3	
7	Std2	
8	Std1	
9		
10		
11		
12		
13		
14		
15		
16		



XI. How to Choose Quantibody® Products?

Species-based arrays:

- <u>Human:</u> QAH-TH-1, QAH-INF-1, QAH-INF-2, QAH-INF-3, QAH-CYT-1, QAH-CYT-2, QAH-MMP-1, QAH-ISO-1, QAH-ANG-1, QAH-ANG-2, QAH-ANG-3, QAH-ANG-1000, QAH-ADI-1, QAH-ADI-2, QAH-CHE-1, QAH-GF-1, QAH-REC-1, QAH-CAA-1000, QAH-CAA-2000, QAH-CAA-3000, QAH-CAA-4000, QAH-CAA-5000, QAH-TH-17
- <u>Mouse:</u> QAM-CYT-1, QAM-CYT-2, QAM-CYT-3, QAM-CYT-4, QAM-CYT-5, QAM-CYT-6, QAM-INF-1, QAM-INT-1, QAM-INT-2, QAM-INT-1000, QAM-CAA-1000, QAM-CYT-Q2000, QAM-CAA-2000, QAM-TH-17
- <u>Rat:</u> QAR-CYT-1, QAR-CYT-2, QAR-CYT-3, QAR-INF-1
- **Porcine:** QAP-CYT-1

Function-based arrays:

- TH1/TH2/TH17 Arrays: QAH-TH-1, QAH-TH-17, QAM-TH-17
- Inflammation Arrays: QAH-INF-1, QAH-INF-2, QAH-INF-3; QAM-INF-1; QAR-INF-1
- Angiogenesis Arrays: QAH-ANG-1, QAH-ANG-2, QAH-ANG-3, QAH-ANG-1000
- MMP Array: QAH-MMP-1
- Immunoglobin Isotype Array: QAH-ISO-1

Cytokine Number-based arrays:

- 240 cytokines: QAH-CAA-5000
- 200 cytokines: QAH-CAA-4000
- 160 cytokines: QAH-CAA-3000
- 120 cytokines: QAH-CAA-2000; QAM-CAA-2000
- 80 cytokines: QAH-CAA-1000; QAM-CAA-1000
- 60 cytokines: QAH-ANG-1000; QAM-CYT-Q2000
- 40 cytokines: QAH-INF-3, QAH-CHE-1, QAH-GF-1, QAH-REC-1, QAM-INF-1, QAM-CYT-4, QAM-CYT-5, QAM-CYT-6, QAH-CYT-4, QAH-CYT-5
- 20-30 cytokines: QAH-ANG-2, QAH-ANG-3, QAM-INT-1000, QAR-CYT-3
- 20 cytokines: QAH-CYT-1, QAM-CYT-1, QAM-CYT-2, QAM-CYT-3, QAM-INT-1
- 10 or less: QAH-TH-1, QAH-INF-1, QAH-INF-2, QAH-ANG-1, QAH-MMP-1, QAH-ADI-1, QAM-INT-2, QAR-CYT-1, QAR-CYT-2, QAR-INF-1, QAH-ISO-1, QAP-CYT-1

Purpose-based array --- Custom Arrays

- Choose from over 400 cytokine pool; Any kind; Any number
- Order slide only or full service in house.

Check our website regularly for updated Quantibody[®] products

Quantibody® Human Immunoglobulin Array 1 20

Note:

Quantibody[®] is the trademark of RayBiotech, Inc.

Cytokine protein arrays are RayBiotech patent-pending technology.

This product is intended for research only and is not to be used for clinical diagnosis. Our produces may not be resold, modified for resale, or used to manufacture commercial products without written approval by RayBiotech, Inc.

Under no circumstances shall RayBiotech be liable for any damages arising out of the use of the materials.

Products are guaranteed for three months from the date of purchase when handled and stored properly. In the event of any defect in quality or merchantability, RayBiotech's liability to buyer for any claim relating to products shall be limited to replacement or refund of the purchase price.

This product is for research use only.



©2010 RayBiotech, Inc.