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Buffer component	Applications in Downstream Bioprocessing
Tris	pH Neutralization, Chromatography
Tris-HCl	Chromatography
Phosphate Buffer	Chromatography, UF/DF
Citrate Buffer	Chromatography, UF/DF
Acetic Acid	pH Acidification, Protein A Chromatography Elution
Citric Acid	pH Acidification, Chromatography, UF/DF
Phosphoric Acid	pH Acidification, Chromatography

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References

1. Bratt, Jennifer, et al. Buffers in Biologics Manufacturing. BioProcess International 15.2(2017).
2. Challener, Cynthia. Behind the Scenes with Buffers. BioPharm International 28.2(2015).
3. Gronemeyer, Petra, et al. Trends in Upstream and Downstream Process Development for Antibody Manufacturing. Bioengineering 1.4(2014).
4. Leveling, Tony. The relationship between pH and Conductivity in a Lithium Contaminated, Deionized Water System. Pbar Note 675: The relationship between pH and Conductivity. Fermilab 2002.
5. West, Jonathan, et al. "PH Measurement." Visual Encyclopedia of Chemical Engineering, University of Michigan College of Engineering, 2018, encyclopedia.che.engin.umich.edu/Pages/ProcessParameters/pHMeasurement/pHMeasurement.html.
6. Manns, Reinhard and Jurgen Schleicher. Information on Conductivity Measurement. JUMO GmbH & Co. KG FAS 624 04.07(2007).
7. 原始文章：[Osmolality as a concentration measurement method for key buffers in bioprocessing](#)