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
Figure 2: SAM265™ methyltransferase assay quantitatively assays Thiopurine S-methyltransferase (EC 2.1.1.67) (TPMT). 0-3µM TPMT was assayed with SAM265, using thiophenol as a substrate. Inset graph shows linear correlation between absorbance and TPMT concentration.

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| EC # | Official Name | Reaction Catalysed |
|--------------|---------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1.16.1.8 | [methionine synthase] reductase | [Methionine synthase]-cob(II) alamin + NADPH + S-adenosyl-L-methionine = [Methionine synthase]-methylcob(II)alamin + AdoHcy + NADP+ |
| 2.1.1.1-157* | S-Adenosylmethionine dependent methyltransferases | S-adenosyl-L-methionine + substrate = AdoHcy + methylated substrate |
| 2.5.1.4 | Adenosylmethionine cyclotransferase | S-adenosyl-L-methionine = MTA + 2-aminobutan-4-olide |
| 2.5.1.16 | Spermidine synthase | S-adenosyl-(5')-3-methylthio-1-propylamine + 1,4-diaminobutane = MTA + spermidine |
| 2.5.1.22 | Spermine synthase | S-adenosylmethioninamine + spermidine = MTA + spermine |
| 2.5.1.23 | Sym-norspermidine synthase | S-adenosylmethioninamine + propane-1,3-diamine = MTA + N-(3-aminopropyl)-1,3-diaminopropane |
| 2.5.1.24 | Discadenine synthase | S-adenosyl-L-methionine + N6-(delta2-isopentenyl)-adenine = MTA + discadenine |
| 2.5.1.25 | tRNA-uridine aminocarboxypropyltransferase | S-adenosyl-L-methionine + tRNA ^{Phe} = MTA + tRNA 3-(3-amino-3-carboxypropyl) uridine |
| 2.5.1.38 | Isonocardicin synthase | S-adenosyl-L-methionine + nocardicin E = MTA + isonocardicin A |
| 2.3.1.161 | Lovastatin nonaketide synthase | acetyl-CoA + malonyl-CoA + NADPH + H ⁺ + S-adenosyl-L-methionine = dihydrononacolin L + CoA + CO ₂ + NADP ⁺ + AdoHcy |
| 3.3.1.2 | Adenosylmethionine hydrolase | S-adenosyl-L-methionine + H ₂ O = MTA + L-homoserine |
| 4.4.1.14 | 1-aminocyclopropane-1-carboxylate synthase | (-)-S-adenosyl-L-methionine = MTA + 1-aminocyclopropane-1-carboxylate |

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

Figure 1: FOCUS™ Glycoprotein isolates a large number of glycoproteins. Jurkat cells were lysed by sonication, centrifuged and the supernatant (CL) loaded onto a FOCUS™ Glycoprotein column. Following incubation the column was centrifuged and the flow through (FT) was collected. The column was washed 5 times (W1-5) and the glycoproteins were sequentially eluted with Glyco Elution Buffer I (E1-2), Glyco Elution Buffer II (E3) and Glyco Elution Buffer III (E4).

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Comparison of crude Jurkat cell lysate (A) and glycoproteins (B) isolated with FOCUS™ Glycoprotein.

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| Cat. # | Description/ size |
|---------|--------------------------------|
| 786-253 | FOCUS™: Glycoprotein/ 10 Preps |