


實驗選擇障礙出現時，不知道該選哪個target好嗎？
經費有限，抗體無法通通買下來時怎麼辦？

iReal現推出 Antibodies Panel，一支抗體的費用，讓您有五支抗體的體驗


EMT Marker



Vimentin N-cadherin

αSMA E-cadherin

Exosome Marker



CD9 ALIX

Panel List




Catalog	Antibodies Panel	Component
IRAP5001	EMT Marker Antibody Panel	Vimentin, Fibronectin, N-cadherin, E-cadherin, Snail Antibody
IRAP5002	ER stress Homeostasis Antibody Panel	IRE1 alpha, CHOP, BiP, PERK, ATF6 Antibody
IRAP5003	Exosome Marker Antibody Panel	CD9, CD63, CD81, HSP70, ALIX Antibody
IRAP5004	ECM Profiling Antibody Panel	Fibronectin, Collagen I, MMP9, MMP1, TIMP1 Antibody
IRAP5005	Myofibroblast / Fibrocyte Antibody Panel	Fibronectin, Collagen I, α-Smooth Muscle Actin, CD34, CD45 Antibody
IRAP5006	Cardiovascular Antibody Panel	CD142 / Tissue Facto, Troponin I, Desmin, α-Smooth Muscle Actin, E-selectin Antibody
IRAP5007	Organelle Marker Antibody Panel	NaK ATPase, EEA1, LAMP2, Lamin B1, PCNA Antibody
IRAP5008	Cell Cycle Marker Antibody Panel	Cyclin A2, Cyclin B1, Cyclin D1, Cyclin E1, p21 Antibody

Western Blot 做得好 實驗沒煩惱

靈敏度提高的最佳利器


BlockPRO Protein-Free Blocking Buffer

- α 無蛋白質：降低 Non Specific Binding
- α 應用廣泛：螢光Western與磷酸化蛋白之首選
- α 提高靈敏度：有效降低背景雜訊
- α 可減少抗體用量



訂購資訊 | BF01-1L | BlockPRO Protein-Free Blocking Buffer | 500 mL*2

LumiFlash™ Ultima長效型冷光呈色劑



1st Ab: Mouse anti-HSP 47 (1:2,000)
2nd Ab: Anti-mouse IgG-HRP (1:20,000)

進階的靈敏度與曝光長度

LumiFlash™ Ultima
冷光試劑組-Ultima 等級偵測靈敏度試驗

K562 cell lysate with 1/2 serial dilution from 50 µg was separated by 12.5% SDS-PAGE and probed by anti HSP 47. All results were exposed to X-ray film for 3 minutes.

冷光CCD照膠系統的最愛


- α 低背景/長時間測試曝光條件最適合
- α Low Picogram靈敏度等級，適合多數實驗室蛋白質表現之應用
- α 減少抗體用量、數據更乾淨

訂購資訊 | LF08-500 | LumiFlash™ Ultima Chemiluminescent Substrate | 1 kit (250 mL solution A + 250 mL solution B)

促銷方案

現購買LF08-500或BF01-1L，即送Protein Marker 250 µL 1支

促銷碼:2022VISDM-009 優惠至2022/11/15止



(貨號:#PM006-0250)



峇祥股份有限公司

THLS-2208-04-022
台北: 02-27851156 新竹: 03-5307592
台中: 04-24710255 高雄: 07-3431735
www.thco.com.tw info@thco.com.tw



LINE@



客戶滿意度



FB



IG