





A Geno Technology, Inc. (USA) brand name

Pre-Diluted BSA Standard

(Cat. # 786-114)



INTRODUCTION

Pre-Diluted BSA Standards are ready to use solutions made in saline buffer containing 0.05% sodium azide. They are provided in six different concentrations ranging from 0.1mg/ml to 1.0mg/ml for running standard calibration in protein estimations.

ITEM(S) SUPPLIED (Cat# 786-114)

Description	Size
Pre-diluted BSA 0.1mg/ml	5ml
Pre-diluted BSA 0.2mg/ml	5ml
Pre-diluted BSA 0.3mg/ml	5ml
Pre-diluted BSA 0.5mg/ml	5ml
Pre-diluted BSA 0.8mg/ml	5ml
Pre-diluted BSA 1.0mg/ml	5ml

STORAGE CONDITIONS

The kit is shipped at ambient temperature. Store refrigerated at 4°C, upon arrival. For long term storage keep at -20°C and is stable for 1 year. Avoid frequent freeze-thaw.


INSTRUCTIONS FOR USE

For use with Non-Interfering Protein Assay kit, use 50µl from each vial for 1.5ml assay format and 5µl per well for 96-well micro-plate assay.

For other types of protein assays use an appropriate volume of the standard as needed.

RELATED PRODUCTS

Download our Protein Assay Handbook



http://info.gbiosciences.com/complete-protein-assay-guide/

For other related products, visit our website at <u>www.GBiosciences.com</u> or contact us.

Last saved: 7/26/2012 CMH

This page is intentionally left blank



www.GBiosciences.com